

RELIABLE | COMMITTED | PROFESSIONAL | INDEPENDENT

ENDOSS



PRODUCT CATALOG 2024

ENDOSCOPY | CLEANING | TRANSPORT


INDEX

PRODUCT INNOVATION	4	Superfine Brush	15
		Superfine Brush.....	15
		ART.NR. 10028.....	15
		Superfine Brush.....	15
		ART.NR. 10023.....	15
		Superfine Brush.....	15
		ART.NR. 10038.....	15
		Push & Pull Brush	16
		Push & Pull Brush.....	16
		ART.NR. 10055.....	16
		Push & Pull Brush.....	16
		ART.NR. 10050.....	16
		Push & Pull Brush.....	16
		ART.NR. 10030.....	16
		Push & Pull Brush.....	16
		ART.NR. 10040.....	17
		Multipoint Brush	17
		Multipoint Brush.....	17
		ART.NR. 10056.....	17
		Multipoint Brush.....	17
		ART.NR. 10058.....	17
		Combi Kit	18
		Channels + Valves Brush Kit.....	18
		ART.NR. 10075.....	18
		Combi Brush	18
		Combi Brush.....	18
		ART.NR. 10065.....	18
		Valves & Ports Brush	19
		2-sided Brush.....	19
		ART.NR. 10060.....	19
		Overview brushes	19
		Custom made brushes	19
ENDOSCOPY	5		
Irrigation tube set	6		
Disposable (24h) irrigation tube set.....	6		
ART.NR. 30081.....	6		
Rinse and insufflation tubing	7		
Disposable (24h) rinse and insufflation tubing.....	7		
ART.NR. 30096.....	7		
Valves for flexible endoscopy	8		
Single use valve set.....	8		
ART.NR. 30077.....	8		
ART.NR. 30078.....	8		
ART.NR. 30075.....	9		
ART.NR. 30074.....	9		
Single use jet adapter Olympus.....	9		
ART.NR. 30055.....	9		
Single use jet adapter Pentax.....	10		
ART.NR. 30054.....	10		
Single use biopsy valve Olympus/Fujinon.....	10		
ART.NR. 30011.....	10		
Single use biopsy valve Pentax.....	10		
ART.NR. 30016.....	10		
Additional products	11		
Single use bite block for adults.....	11		
ART.NR. 30025.....	11		
Reusable metal adapter UCR.....	11		
ART.NR. 30060.....	11		
Examination pants colonoscopy.....	12		
ART.NR. 55015.....	12		
Water-resistant gown.....	12		
ART.NR. 55012.....	12		
CLEANING	13		
Single en Double Brush	14		
Single Brush.....	14		
ART.NR. 10010.....	14		
Double Brush.....	14		
ART.NR. 10020.....	14		


TRANSPORT	20
Tip protection	21
Mesh Protector small red.....	21
ART.NR. 22026.....	21
Mesh Protector large blue.....	21
ART.NR. 22027.....	21
Scope Guard small.....	22
ART.NR. 22006.....	22
Scope Guard medium.....	22
ART.NR. 22007.....	22
Scope Guard large.....	22
ART.NR. 22008.....	22
Smiley covers	23
Smiley cover for transport tray 60x40.....	23
ART.NR. 22014.....	23
Smiley cover for transport tray 52x42.....	23
ART.NR. 22017.....	23
Smiley cover for transport tray 40x30.....	23
ART.NR. 22013.....	23
Smiley cover for transport tray 51x51.....	24
ART.NR. 22011.....	24
Smiley cover for transport tray 48x25x10 (DIN).24	
ART.NR. 22019.....	24
ABOUT ENDOSS B.V.	25
OUR SALES PARTNERS	25

PRODUCT INNOVATION

Designing and developing products is an important part of ENDOSS. By listening carefully to the users of our products, we strive to continuously improve them. Equally important are the new ideas that emerge from conversations we and our distributors have in hospitals and clinics.



2Catch Polyp Trap
in development



Tip Protector
in development



ENDOSCOPY

Endoss develops, designs, and distributes medical disposables. We are committed to improve patient safety and infection prevention. We do this by offering products for CSSD and Endoscopy, our fields of expertise and focus.

Patient safety: that is what our goal is. We strive to improve the safety of the patient and take it to the next level. For this reason most of our disposable products are sterilised en packed separately in user friendly packings. All our products are labeled in accordance with the MDR, or other applicable standards like FDA. Each packing is provided with a GS1 barcode or datamatrix (UDI).

IRRIGATION TUBE SET

Endoss has disposable irrigation tube sets available for different types of brands flushing pumps. In order to prevent micro-organism to get in contact with the sterile water, a filter is placed in the cap. Research has shown that after 24 hour the water is still bacteria free.

Other advantages of this product are:

- Disposable (24 hour)
- Universal cap
- Suitable for the most familiar flushing pumps
- Filter in the cap, sterile water remains sterile
- Back-flow valve in the connector




Disposable (24h) irrigation tube set
For Erbe and Olympus flushing pumps

ART.NR. 30081
VE 4x10

RINSE AND INSUFFLATION TUBING

Besides irrigation tube sets Endoss has disposable rinse and insufflation tube sets. Suitable for the most well-known brands endoscopes and CO2 insufflators. With an long extension for the CO2 it is no longer necessary to purchase a separate CO2 tube.



Other advantages of this product:

- Disposable (24 hour)
- Universal cap
- Suitable for the most well-known brands endoscopes and CO2 insufflators
- Filter in the cap, sterile water remains sterile
- Back-flow valve in the CO2 tube
- Standard for Olympus, systems for other brands available upon request



Disposable (24h) rinse and insufflation tubing
 For Olympus with filter in the cap and long CO2 extension

ART.NR. 30096
 VE 4x10

VALVES FOR FLEXIBLE ENDOSCOPY

Reusable valves for endoscopes require a difficult and precise cleaning procedure that needs to be followed according to the manufacturer’s instructions. With our single use valves manually cleaning becomes unnecessary. The single use valves contribute to the patient safety and ensure an optimized infection prevention.

Features:

- Available as sets but also as separate products
- Less effort for the CSSD employers
- Increase patient safety
- Available for all well-known brands endoscopes



Single use valve set

With suction, air/water, biopsy Valve & auxiliary water connector
Suitable for Olympus endoscopes

ART.NR. 30077
VE 1x50



Single use valve set

Suction, air/water, biopsy Valve & auxiliary water connector
Suitable for Pentax endoscopes

ART.NR. 30078
VE 1x50



Single use valve set

Suction, air/water, biopsy Valve & auxiliary water connector for Olympus endoscopes

ART.NR. 30075
VE 1x50



Single use valve set air/water

For Olympus endoscopes

ART.NR. 30074
VE 1x50



Single use auxiliary adapter

For Olympus endoscopes

ART.NR. 30055
VE 1x50

| Endoscopy - Valves for flexible endoscopy



Single use auxiliary adapter
For Pentax endoscopes

ART.NR. 30054
VE 1x50



Single use biopsy valve
For Olympus/Fujinon endoscopes

ART.NR. 30012
VE 4x100



Single use biopsy valve
For Pentax endoscopes

ART.NR. 30016
VE 4x100

ADDITIONAL PRODUCTS

Besides disposable valves, irrigation tubes and rinse and insufflation tubes we offer various other products that improve patient safety or facilitates the process.

These include:

- Bite blocks, for both adults and pediatric examinations
- Metal UCR adapter
- Metal bottle holder
- Polyp trap (under development)



Single use bite block for adults
With patient friendly elastic headband

ART.NR. 30025
VE 4x100



Reusable metal adapter UCR
Enables to attach the disposable rinse and insufflation tube on the UCR insufflator

ART.NR. 30060
VE 1

| **Endoscopy - Additional products**

During an examination it is necessary that patient and staff feel secure and safe. With our clothing we offer safety and discretion for the patient and staff/doctors.

Features:

- Single use clothing
- Different sizes available
- Comfortable



Examination pants colonoscopy

With an opening at the back with elastic band

ART.NR. 55015
VE 1x100

The staff needs to be able to work in a safe way without being exposed to contaminated fluids and physical material.

Features:

- Work safely
- Protection during the examination
- Protection during the precleaning
- Made of high quality material



Water-resistant gown

With long sleeves and practical thumb openings

ART.NR. 55012
VE 1x100



CLEANING

Endoss offers various solutions for a solid pre cleaning of flexible endoscopes. We offer brushes for the pre cleaning of the channels of flexible endoscopes.

Over the 10 years Endoss has sold their brushes throughout Europe. In this period we produced and delivered millions of brushes. We have a wide range of brushes for all types and brands of endoscopes.


SINGLE AND DOUBLE BRUSH

Various methods and endoscopes require different types of brushes. The Single and the Double Brush are brushes that can be used for the pre cleaning of the channels of the endoscopes. The Single Brush can be pushed in and pulled back. The Double Brush can be inserted in the channel of the endoscope and pulled through the whole channel.

Brushes can be custom made if your desired type is not in our product range.


Features:

- 5mm diameter brush
- User friendly
- Suitable for pull through and pulling back the brush



Single Brush L230cm

ART.NR. 10010
VE 4x100



Double Brush L230cm


ART.NR. 10020
VE 4x100

SUPERFINE BRUSH

Very fine brushes make it possible to clean the smallest channels.


Features:

- Suitable for channel diameters between 1.0 and 1.5mm
- Individually packaged and equipped with GS1 code





Superfine Brush L110cm

ART.NR. 10028
VE 1x50



Superfine Brush L230cm

ART.NR. 10023
VE 1x50



Superfine Brush L380cm

ART.NR. 10038
VE 1x50

PUSH & PULL BRUSH

The Push & Pull brushes of Endoss are brushes that offer an excellent quality and cleaning effectiveness of LOG4. The brushes need to be pushed into the channel and pulled through.


Features:

- Significantly better cleaning results LOG4
- Available in various diameters and lengths
- One usage of the Push & Pull is three times better than traditionally pre cleaning
- Packed per piece and marked with a GS1 code



Push&Pull Brush L230cm

ART.NR. 10055
VE 4x100



Push&Pull Brush L160cm



ART.NR. 10050
VE 4x100



Push&Pull Brush L160cm

ART.NR. 10030
VE 4x100

| **Cleaning - Push & Pull and Multipoint brush**



Push&Pull Brush L120cm



ART.NR. 10040
VE 4x100

MULTIPOINT BRUSH

Our Multipoint brush exists of three brushes each 20mm long and 6.0mm in diameter. This guarantees an improved cleaning efficacy



Features:

- Three brushes
- Easy in use
- Length 230cm



Multipoint Brush L230cm

ART.NR. 10056
VE 4x100



Multipoint Brush L230cm

ART.NR. 10058
VE 4x100

COMBI KIT

A lot of people use our Push & Pull brush and the 2-sided brush that is why we offer a kit where both products are included.

Advantages:

- Saving space on the department (less boxes)
- Less article numbers
- Less risk of reusing the brushes



Channels + Valves Brush Kit

10055 - Push&Pull Brush L230cm 5.0mm
 10060 - 2-sided brush 5mm + 11mm


ART.NR. 10075
 VE 4x150

COMBI BRUSH

The Combi Brush is a convenient channel brush with a valve brush at the end.

Features:

- Handle for extra grip during retrieval
- Channel and valve brush in one



Combi Brush 6mm+11mm

ART.NR. 10065
 VE 4x100

VALVES & PORTS BRUSH

The 2-sided brush is intended for cleaning the valves and ports of an endoscope.

Features:

- Handle for additional grip
- Packaged together in a dispenser box
- Two brush diameters in one



2-Sided Brush 5mm+11mm

ART.NR. 10060
VE 4x100

OVERVIEW BRUSHES

WIRE LENGTH	Ø BRUSH 1.5 MM	Ø BRUSH 3.0 MM	Ø BRUSH 4.0 MM	Ø BRUSH 5.0 MM	Ø BRUSH 6.0 MM	Ø BRUSH 11.0 MM
110 CM	• SUPERFINE 10028					
120 CM		• PUSH&PULL 10040				
160 CM		• PUSH&PULL 10030		• PUSH&PULL 10050		
230 CM	• SUPERFINE 10023		• MULTIPOINT 10056	• SINGLE 10010 • DOUBLE 10020 • PUSH&PULL 10055	• MULTIPOINT 10058 • COMBI 10065	
380 CM	• SUPERFINE 10038					
VALVES				• 2-SIDED 10060		• 2-SIDED 10060 • COMBI 10065

CUSTOM MADE BRUSHES

Are you looking for a specific brush, but it is not in our product range yet? No problem, we also offer brushes that are custom made.



TRANSPORT


Endoss provides a range of solutions for clearly marking both dirty and clean transport processes. Additionally, we offer products designed to safeguard endoscopes during transportation.

TIP PROTECTION

Tip protectors ensure that endoscopes or medical devices can be protected during transport and handling. We offer different types in our product range from various type of materials.

Features:

- Sterilized
- Suitable for all sizes of endoscopes
- A moisture permeable type
- A type that can be used as a sponge as well



Mesh Protector small red

ART.NR. 22026
VE 4x100



Mesh Protector large blue

ART.NR. 22027
VE 4x100



Scope Guard small

ART.NR. 22006
VE 4x50



Scope Guard medium

ART.NR. 22007
VE 4x50



Scope Guard Large
Also suitable for cleaning the outside of the endoscope

ART.NR. 22008
VE 4x50

SMILEY COVERS

For a clear marking of dirty and clean transport we offer smiley covers. We have all types of sizes available for various transport trays.

Advantages:

- Sterilized
- Available in different sizes
- Smiley on both sides



Smiley cover for transport tray 60x40

ART.NR. 22014
VE 4x50



Smiley cover for transport tray 52x42

ART.NR. 22017
VE 4x50



Smiley cover for transport tray 40x30

ART.NR. 22013
VE 2x50



Smiley cover for transport tray 51x51

ART.NR. 22011
VE 1x300



Smiley cover for transport tray 48x25x10 (DIN)

ART.NR. 22019
VE 4x50


ABOUT ENDOSS B.V.

ENDOSS is a leading Dutch manufacturer and distributor of medical devices, specializing in flexible endoscopy, patient safety, and infection prevention. Since 2008, we have been supplying our high-quality products to hospitals, clinics, and practices worldwide.

We operate our own branches in the Netherlands and Germany, while in other countries, we work with exclusive distributors.

ENDOSS stands for reliability, safety, user-friendliness, and innovation.

HERE YOU FIND OUR GLOBAL SALES PARTNERS



THE NETHERLANDS

ENDOSS BV
Minervum 7166b
4817 ZN Breda, The Netherlands

+31 (0)76 53 09 268
info@endoss.nl
www.endoss.nl

GERMANY

ENDOSS GmbH
Viktoriastraße 10
46483 Wesel, Germany

+49 (0)281 97 45 3001
office@endoss.de
www.endoss.de