

RELIABLE | COMMITTED | PROFESSIONAL | INDEPENDENT

# ENDOSS



## PRODUCT CATALOG 2025

ENDOSCOPY | CLEANING | TRANSPORT

# INDEX

<b>PRODUCT INNOVATION</b>	<b>4</b>	<b>Superfine Brush</b>	<b>15</b>
		Superfine Brush.....	15
		ART.NR. 10028.....	15
		Superfine Brush.....	15
		ART.NR. 10023.....	15
		Superfine Brush.....	15
		ART.NR. 10038.....	15
		<b>Push &amp; Pull Brush</b>	<b>16</b>
		Push & Pull Brush.....	16
		ART.NR. 10055.....	16
		Push & Pull Brush.....	16
		ART.NR. 10050.....	16
		Push & Pull Brush.....	16
		ART.NR. 10030.....	16
		Push & Pull Brush.....	16
		ART.NR. 10040.....	17
		<b>Multipoint Brush</b>	<b>17</b>
		Multipoint Brush.....	17
		ART.NR. 10056.....	17
		Multipoint Brush.....	17
		ART.NR. 10058.....	17
		<b>Combi Kit</b>	<b>18</b>
		Channels + Valves Brush Kit.....	18
		ART.NR. 10075.....	18
		<b>Combi Brush</b>	<b>18</b>
		Combi Brush.....	18
		ART.NR. 10065.....	18
		<b>Valves &amp; Ports Brush</b>	<b>19</b>
		2-sided Brush.....	19
		ART.NR. 10060.....	19
		<b>Overview brushes</b>	<b>19</b>
		<b>Custom made brushes</b>	<b>19</b>
<b>ENDOSCOPY</b>	<b>5</b>		
<b>Irrigation tube set</b>	<b>6</b>		
Disposable (24h) irrigation tube set.....	6		
ART.NR. 30081.....	6		
<b>Rinse and insufflation tubing</b>	<b>7</b>		
Disposable (24h) rinse and insufflation tubing.....	7		
ART.NR. 30096.....	7		
<b>Valves for flexible endoscopy</b>	<b>8</b>		
Single use valve set.....	8		
ART.NR. 30077.....	8		
ART.NR. 30078.....	8		
ART.NR. 30075.....	9		
ART.NR. 30074.....	9		
Single use jet adapter Olympus.....	9		
ART.NR. 30055.....	9		
Single use jet adapter Pentax.....	10		
ART.NR. 30054.....	10		
Single use biopsy valve Olympus/Fujinon.....	10		
ART.NR. 30011.....	10		
Single use biopsy valve Pentax.....	10		
ART.NR. 30016.....	10		
<b>Additional products</b>	<b>11</b>		
Single use bite block for adults.....	11		
ART.NR. 30025.....	11		
Reusable metal adapter UCR.....	11		
ART.NR. 30060.....	11		
Examination pants colonoscopy.....	12		
ART.NR. 55015.....	12		
<b>CLEANING</b>	<b>13</b>		
<b>Single en Double Brush</b>	<b>14</b>		
Single Brush.....	14		
ART.NR. 10010.....	14		
Double Brush.....	14		
ART.NR. 10020.....	14		


<b>TRANSPORT</b>	<b>20</b>
<b>Tip protection</b>	<b>21</b>
Mesh Protector small red.....	21
ART.NR. 22026.....	21
Mesh Protector large blue.....	21
ART.NR. 22027.....	21
Scope Guard small.....	22
ART.NR. 22006.....	22
Scope Guard medium.....	22
ART.NR. 22007.....	22
Scope Guard large.....	22
ART.NR. 22008.....	22
<b>Smiley covers</b>	<b>23</b>
Smiley cover for transport tray 60x40.....	23
ART.NR. 22014.....	23
Smiley cover for transport tray 52x42.....	23
ART.NR. 22017.....	23
Smiley cover for transport tray 40x30.....	23
ART.NR. 22013.....	23
Smiley cover for transport tray 51x51.....	24
ART.NR. 22011.....	24
Smiley cover for transport tray 48x25x10 (DIN).24	
ART.NR. 22019.....	24
<b>ABOUT ENDOSS B.V.</b>	<b>25</b>
<b>OUR SALES PARTNERS</b>	<b>25</b>

**PRODUCT INNOVATION**

---


Designing and developing products is an important part of ENDOSS. By listening carefully to the users of our products, we strive to continuously improve them. Equally important are the new ideas that emerge from conversations we and our distributors have in hospitals and clinics.

---



**2Catch Polyp Trap**  
in development

---



**Tip Protector**  
in development

---



# ENDOSCOPY

Endoss develops, designs, and distributes medical disposables. We are committed to improve patient safety and infection prevention. We do this by offering products for CSSD and Endoscopy, our fields of expertise and focus.

Patient safety: that is what our goal is. We strive to improve the safety of the patient and take it to the next level. For this reason most of our disposable products are sterilised en packed separately in user friendly packings. All our products are labeled in accordance with the MDR, or other applicable standards like FDA. Each packing is provided with a GS1 barcode or datamatrix (UDI).

## IRRIGATION TUBE SET

Endoss has disposable irrigation tube sets available for different types of brands flushing pumps. In order to prevent micro-organism to get in contact with the sterile water, a filter is placed in the cap. Research has shown that after 24 hour the water is still bacteria free.

Other advantages of this product are:

- Disposable (24 hour)
- Universal cap
- Suitable for the most familiar flushing pumps
- Filter in the cap, sterile water remains sterile
- Back-flow valve in the connector




**Disposable (24h) irrigation tube set**  
For Erbe and Olympus flushing pumps

ART.NR. 30081  
VE 4x10

**RINSE AND INSUFFLATION TUBING**

Besides irrigation tube sets Endoss has disposable rinse and insufflation tube sets. Suitable for the most well-known brands endoscopes and CO2 insufflators. With an long extension for the CO2 it is no longer necessary to purchase a separate CO2 tube.



Other advantages of this product:

- Disposable (24 hour)
- Universal cap
- Suitable for the most well-known brands endoscopes and CO2 insufflators
- Filter in the cap, sterile water remains sterile
- Back-flow valve in the CO2 tube
- Standard for Olympus, systems for other brands available upon request



**Disposable (24h) rinse and insufflation tubing**  
 For Olympus with filter in the cap and long CO2 extension

ART.NR. 30096  
 VE 4x10

**VALVES FOR FLEXIBLE ENDOSCOPY**

Reusable valves for endoscopes require a difficult and precise cleaning procedure that needs to be followed according to the manufacturer’s instructions. With our single use valves manually cleaning becomes unnecessary. The single use valves contribute to the patient safety and ensure an optimized infection prevention.

Features:

- Available as sets but also as separate products
- Less effort for the CSSD employers
- Increase patient safety
- Available for all well-known brands endoscopes



**Single use valve set**

With suction, air/water, biopsy Valve & auxiliary water connector  
Suitable for Olympus endoscopes

ART.NR. 30077  
VE 1x50



**Single use valve set**

Suction, air/water, biopsy Valve & auxiliary water connector  
Suitable for Pentax endoscopes

ART.NR. 30078  
VE 1x50



**Single use valve set**

Suction, air/water, biopsy Valve & auxiliary water connector for Olympus endoscopes

ART.NR. 30075  
VE 1x50



**Single use valve set air/water**

For Olympus endoscopes

ART.NR. 30074  
VE 1x50



**Single use auxiliary adapter**

For Olympus endoscopes

ART.NR. 30055  
VE 1x50

| Endoscopy - Valves for flexible endoscopy



**Single use auxiliary adapter**  
For Pentax endoscopes

ART.NR. 30054  
VE 1x50



**Single use biopsy valve**  
For Olympus/Fujinon endoscopes

ART.NR. 30012  
VE 4x100



**Single use biopsy valve**  
For Pentax endoscopes

ART.NR. 30016  
VE 4x100

## ADDITIONAL PRODUCTS

Besides disposable valves, irrigation tubes and rinse and insufflation tubes we offer various other products that improve patient safety or facilitates the process.

These include:

- Bite blocks, for both adults and pediatric examinations
- Metal UCR adapter
- Metal bottle holder
- Polyp trap (under development)



**Single use bite block for adults**  
With patient friendly elastic headband

ART.NR. 30025  
VE 4x100



**Reusable metal adapter UCR**  
Enables to attach the disposable rinse and insufflation tube on the UCR insufflator

ART.NR. 30060  
VE 1

## | Endoscopy - Additional products

During an examination it is necessary that patient and staff feel secure and safe. With our clothing we offer safety and discretion for the patient and staff/doctors.

Features:

- Single use clothing
- Different sizes available
- Comfortable



**Examination pants colonoscopy**  
With an opening at the back with elastic band

ART.NR. 55015  
VE 1x100



## CLEANING

Endoss offers various solutions for a solid pre cleaning of flexible endoscopes. We offer brushes for the pre cleaning of the channels of flexible endoscopes.

Over the 10 years Endoss has sold their brushes throughout Europe. In this period we produced and delivered millions of brushes. We have a wide range of brushes for all types and brands of endoscopes.


**SINGLE AND DOUBLE BRUSH**

Various methods and endoscopes require different types of brushes. The Single and the Double Brush are brushes that can be used for the pre cleaning of the channels of the endoscopes. The Single Brush can be pushed in and pulled back. The Double Brush can be inserted in the channel of the endoscope and pulled through the whole channel.

Brushes can be custom made if your desired type is not in our product range.


Features:

- 5mm diameter brush
- User friendly
- Suitable for pull through and pulling back the brush



**Single Brush L230cm**

ART.NR. 10010  
VE 4x100



**Double Brush L230cm**


ART.NR. 10020  
VE 4x100

**SUPERFINE BRUSH**

Very fine brushes make it possible to clean the smallest channels.


Features:

- Suitable for channel diameters between 1.0 and 1.5mm
- Individually packaged and equipped with GS1 code





**Superfine Brush L110cm**

ART.NR. 10028  
VE 1x50



**Superfine Brush L230cm**

ART.NR. 10023  
VE 1x50



**Superfine Brush L380cm**

ART.NR. 10038  
VE 1x50

**PUSH & PULL BRUSH**

The Push & Pull brushes of Endoss are brushes that offer an excellent quality and cleaning effectiveness of LOG4. The brushes need to be pushed into the channel and pulled through.


Features:

- Significantly better cleaning results LOG4
- Available in various diameters and lengths
- One usage of the Push & Pull is three times better than traditionally pre cleaning
- Packed per piece and marked with a GS1 code



**Push&Pull Brush L230cm**

ART.NR. 10055  
VE 4x100



**Push&Pull Brush L160cm**

ART.NR. 10050  
VE 4x100



**Push&Pull Brush L160cm**

ART.NR. 10030  
VE 4x100

| **Cleaning - Push & Pull and Multipoint brush**



**Push&Pull Brush L120cm**

ART.NR. 10040  
VE 4x100

**MULTIPOINT BRUSH**

Our Multipoint brush exists of three brushes each 20mm long and 6.0mm in diameter. This guarantees an improved cleaning efficacy

Features:

- Three brushes
- Easy in use
- Length 230cm



**Multipoint Brush L230cm**

ART.NR. 10056  
VE 4x100



**Multipoint Brush L230cm**

ART.NR. 10058  
VE 4x100

**COMBI KIT**

A lot of people use our Push & Pull brush and the 2-sided brush that is why we offer a kit where both products are included.

Advantages:

- Saving space on the department (less boxes)
- Less article numbers
- Less risk of reusing the brushes



**Channels + Valves Brush Kit**

10055 - Push&Pull Brush L230cm 5.0mm  
 10060 - 2-sided brush 5mm + 11mm


ART.NR. 10075  
 VE 4x150

**COMBI BRUSH**

The Combi Brush is a convenient channel brush with a valve brush at the end.

Features:

- Handle for extra grip during retrieval
- Channel and valve brush in one



**Combi Brush 6mm+11mm**

ART.NR. 10065  
 VE 4x100

**VALVES & PORTS BRUSH**

The 2-sided brush is intended for cleaning the valves and ports of an endoscope.

Features:

- Handle for additional grip
- Packaged together in a dispenser box
- Two brush diameters in one



**2-Sided Brush 5mm+11mm**

ART.NR. 10060  
VE 4x100

**OVERVIEW BRUSHES**

WIRE LENGTH	Ø BRUSH 1.5 MM	Ø BRUSH 3.0 MM	Ø BRUSH 4.0 MM	Ø BRUSH 5.0 MM	Ø BRUSH 6.0 MM	Ø BRUSH 11.0 MM
110 CM	• SUPERFINE 10028					
120 CM		• PUSH&PULL 10040				
160 CM		• PUSH&PULL 10030		• PUSH&PULL 10050		
230 CM	• SUPERFINE 10023		• MULTIPOINT 10056	• SINGLE 10010 • DOUBLE 10020 • PUSH&PULL 10055	• MULTIPOINT 10058 • COMBI 10065	
380 CM	• SUPERFINE 10038					
VALVES				• 2-SIDED 10060		• 2-SIDED 10060 • COMBI 10065

**CUSTOM MADE BRUSHES**

Are you looking for a specific brush, but it is not in our product range yet? No problem, we also offer brushes that are custom made.



# TRANSPORT


Endoss provides a range of solutions for clearly marking both dirty and clean transport processes. Additionally, we offer products designed to safeguard endoscopes during transportation.

**TIP PROTECTION**

Tip protectors ensure that endoscopes or medical devices can be protected during transport and handling. We offer different types in our product range from various type of materials.

Features:

- Sterilized
- Suitable for all sizes of endoscopes
- A moisture permeable type
- A type that can be used as a sponge as well



**Mesh Protector small red**

ART.NR. 22026  
VE 4x100



**Mesh Protector large blue**

ART.NR. 22027  
VE 4x100



**Scope Guard small**

ART.NR. 22006  
VE 4x50



**Scope Guard medium**

ART.NR. 22007  
VE 4x50



**Scope Guard Large**  
Also suitable for cleaning the outside of the endoscope

ART.NR. 22008  
VE 4x50

**SMILEY COVERS**

For a clear marking of dirty and clean transport we offer smiley covers. We have all types of sizes available for various transport trays.

Advantages:

- Sterilized
- Available in different sizes
- Smiley on both sides



**Smiley cover for transport tray 60x40**

ART.NR. 22014  
VE 4x50



**Smiley cover for transport tray 52x42**

ART.NR. 22017  
VE 4x50



**Smiley cover for transport tray 40x30**

ART.NR. 22013  
VE 2x50



Smiley cover for transport tray 51x51

ART.NR. 22011  
VE 1x300



Smiley cover for transport tray 48x25x10 (DIN)

ART.NR. 22019  
VE 4x50

## ABOUT ENDOSS B.V.

---


ENDOSS is a leading Dutch manufacturer and distributor of medical devices, specializing in flexible endoscopy, patient safety, and infection prevention. Since 2008, we have been supplying our high-quality products to hospitals, clinics, and practices worldwide.

We operate our own branches in the Netherlands and Germany, while in other countries, we work with exclusive distributors.

ENDOSS stands for reliability, safety, user-friendliness, and innovation.

## HERE YOU FIND OUR GLOBAL SALES PARTNERS

---



### THE NETHERLANDS

ENDOSS B.V.  
Ambachtsweg 4  
4854 MK Bavel  
The Netherlands  
[www.endoss.nl](http://www.endoss.nl)

**Customer service**  
**+31-76-303 0343**  
**[order@endoss.nl](mailto:order@endoss.nl)**

General questions  
+31-76-530 9268  
[info@endoss.nl](mailto:info@endoss.nl)